


Experiment settings

Mask

Peak picking

Geometry fitting

Cake &amp; integration



X: -127.3504

Y: 763.4291

Value: n/a

How to

The target is to identify at least 2 rings by location and number. Then to extract all peaks automatically.

Click on the ring you want to select. Usually it is the first one, else update it's number in the list of the picked rings.

Use the recalibration tool to extract more peaks automatically.

Pick peaking

Mode:



Ring



Single pick

Picked rings

	Name	Peaks	Ring number	
	a	227	1	
	b	177	2	
	c	146	3	
	d	132	4	

Undo extract rings

Redo

Recalibrate

Max rings to extract:

4

Number of peak per degree:

1.00

Extract

Next &gt;